

Product Datasheet

ADRB2 antibody (orb19116)

Description	Goat polyclonal antibody to ADRB2
Species/Host	Goat
Reactivity	Bovine, Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ADRB2R / ADRB2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 48kDa band observed in human liver lysates (calculated MW of 46.6kDa according to NP_000015.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	HQGTVPDNDISQ
Clonality	Polyclonal
MW	46.6
Entrez	154
Dilution Range	ELISA: 1:16000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb19116
(0.2ug/ml)
staining of
human li...