


Product Datasheet

EPB41L3 antibody (orb19117)

Description	Goat polyclonal antibody to EPB41L3
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	EPB41L3 / DAL1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Approx 150kDa band observed in human brain lysates (calculated MW of 121kDa according to NP_036439.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	HKETEITPEDGED
Clonality	Polyclonal
MW	121
Entrez	23136
Dilution Range	ELISA: 1:2000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb19117
(1µg/ml)
staining of
human brai...



Western blot
analysis of
Human
Cerebellu...