


Product Datasheet

SLC9A3R2 antibody (orb19121)

Description	Goat polyclonal antibody to SLC9A3R2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	NHERF2 (isoform a)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: In paraffin embedded Human Kidney shows strong membranous staining in glomerulus-. Recommended concentration, 1-2ug/ml.WB: Approx 38kDa band observed in human kidney lysates (calculated MW of 37.4kDa according to BAA33216.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	DGSAWKQDPFQ
Clonality	Polyclonal
MW	37.4
Entrez	9351
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml, IHC-P: 1-2µg/ml
Expiration Date	12 months from date of receipt.



orb19121
(0.3ug/ml)
staining of
human ki...



orb19121
(3ug/ml)
staining of
paraffin e...