



Product Datasheet

GATA1 antibody (orb19125)


Description	Goat polyclonal antibody to GATA1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Target	GATA1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 48kDa band observed in lysates of human liver, placenta and peripheral blood mononucleocytes (calculated MW of 42.8kDa according to NP_002040.1. Recommended concentration: 0.3-1.0 µg/ml.
Protein Sequence	DAEAYRHSPVFQ
Clonality	Polyclonal
MW	42.8
Entrez	2623
Dilution Range	ELISA: 1:16000, WB: 0.3-1 µg/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of paraformalde...



Immunofluorescence analysis of paraforma...



Immunofluorescence analysis of paraforma...