


Product Datasheet

FBXO43 antibody (orb19137)

Description	Goat polyclonal antibody to FBXO43
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	EMI2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: This product has been successfully used on Rat Pachytene Spermatocytes showing a 71kDa band (Perrard et al, PLoS One. 2009 Oct 5;4(10):e7237. PMID: 19802389). Recommended concentration: 1-3 µg/ml.Experiment Notes: This product has been successfully used in Immunofluorescence and Flow Cytometry on Rat meiotic Testis cells: (Perrard et al, PLoS One. 2009 Oct 5;4(10):e7237. PMID: 19802389). Recommended concentration: 5-10ug/ml.
Protein Sequence	PGSAQSKRNLKRL
Clonality	Polyclonal
MW	78.4
Entrez	286151
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Ovary lys...