

Product Datasheet

VPS25 Antibody (orb19138)

Description	Goat polyclonal antibody to VPS25	Image not found or type unknown
Reactivity	Canine, Human, Mouse, Rat	orb19138 (0.5ug/ml) staining of human li...
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	Image not found or type unknown
Target	VPS25	orb19138 (10ug/ml) staining of paraffin ...
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
Note	For research use only	
Application notes	WB: Approx. 22kDa band observed in Rat Kidney lysates (calculated MW of 20.7kDa according to human NP_115729.1). Recommended concentration: 0.3-1µg/ml. Primary incubation was 1 hour.. IHC: Paraffin embedded Human Kidney. Recommended concentration: 6-8µg/ml.	
Protein Sequence	QPNVDTRQKQ	
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml	
Expiration Date	12 months from date of receipt.	