



Product Datasheet

SMYD5 antibody (orb191485)


Description	Rabbit polyclonal antibody to SMYD5
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human SMYDSMYD5
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48kDa
Uniprot ID	Q6GMV2
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100, IF/ICC: 1:50-100, IP: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemical analysis of formalin...



Immunofluorescence analysis of A549 cell...