



Product Datasheet

ROR2 antibody (orb191497)


Description	ROR2 antibody
Species/Host	Rabbit
Reactivity	Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	Recombinant protein
Target	ROR2
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Uniprot ID	Q01974
NCBI	004551 , 004560
Entrez	4920
Dilution Range	WB: 1:500, IHC-P: 1:200, IF/ICC: 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of breast cancer t...



Immunofluorescence analysis of pig kidney...



IHC-P image of pig kidney tissue using R...