

Product Datasheet

CXCL11 antibody (orb191499)

Description

CXCL11 antibody

Species/Host

Rabbit

Reactivity

Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, ICC, IF, IHC-P, WB

Immunogen

KLH conjugated synthetic peptide derived from human CXCL11. Please contact us for the exact immunogen sequence. The peptide is available as orb393929.

Target

CXCL11

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

10 kDa

Uniprot ID

[Q8R392](#), [Q7TNL0](#)

NCBI

[005409](#), [005400](#), [41](#), [11](#)

Entrez

[6373](#)


Dilution Range

WB: 2µg/ml, IHC-P:1:200, IF/ICC:1:200


Expiration Date

12 months from date of receipt.


Q8 R392 Q7TNL0 Merge



IF image of rat brain tissue using anti-...



Immunofluorescence image of rat brain ti...



Immunofluorescence analysis of rat brain...