

Product Datasheet

VEGF antibody (orb191500)

Description

VEGF antibody


Species/Host

Rabbit

Reactivity


Bovine, Canine, Gallus, Guinea pig, Human, Mouse, Porcine, Rat, Sheep

Conjugation

Unconjugated

Tested
Applications

ELISA, ICC, IF, IHC-P, WB



Western blot analysis of rat brain tissue...

Immunogen

KLH conjugated synthetic peptide derived from human VEGF. Please contact us for the exact immunogen sequence. The peptide is available as orb374829.

Target


VEGF

Form/Appearance

10 mM PBS, 0.02% sodium azide, 10mg/ml rAlbumin, 25% glycerol

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml



IHC-P image of mouse skin tissue using a...

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles

Note

For research use only

Isotype


IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography



Immunofluorescence analysis of mouse skin...

Uniprot ID
[P15692](#), [Q00731](#), [P16612](#)
NCBI
[001020537](#), [001025366](#)
Entrez
[7422](#)
Dilution Range

IF/ICC:1:100, IHC-P: 1:50, WB: 1:200

Expiration Date

12 months from date of receipt.