

Product Datasheet

IFT88 antibody (orb19151)

Description

Goat polyclonal antibody to IFT88

Species/Host

Goat

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, IF, WB

Applications
Target

IFT88 / Polaris

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:32000. IHC: In paraffin embedded Human Kidney shows staining of the brush border in PCT. Recommended concentration, 3ug/ml. WB: Preliminary experiments in Human Kidney, Liver and Lung lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates

Protein Sequence

KKRIDEEDFADEE

Clonality

Polyclonal

MW


94.3; 93.2

Entrez
8100
Dilution Range

ELISA: 1:128000, WB: 1-3ug/ml, IHC-P: 10ug/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of Human Ki...