



Product Datasheet

ALDH1A1 antibody (orb19153)

Description	Goat polyclonal antibody to ALDH1A1
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ALDH1A1 (Internal)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 55kDa band observed in Human Liver lysates (calculated MW of 54.9kDa according to NP_000680.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	RELGEYGFHEYTE
Clonality	Polyclonal
MW	54.9
Entrez	216
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Liver Lys...



Immunohistochemical staining of Human Te...



Western blot analysis of Human Liver Lys...