



Product Datasheet

CDKN2A antibody (orb19161)

Description	Goat polyclonal antibody against CDKN2A
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	CDKN2A (isoform 1) / p16INK4a
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx. 16-18kDa band observed in lysates of the cell lines A431 and HeLa (calculated MW of 16.5kDa according to NP_000068.1). Recommended concentration: 0.1-1 µg/ml.
Protein Sequence	HARIDAAEGPSDIPD
Clonality	Polyclonal
MW	16.5; 17.9
Entrez	1029
Dilution Range	ELISA: 1:32000, WB: 0.1-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb19161
(1µg/ml)
staining of
HeLa Lysat...



orb19161
(5µg/ml)
staining of
paraffin ...