

Product Datasheet

CRMP-2 (Phospho-Thr509) Antibody (orb191641)

Description	CRMP-2 (Phospho-Thr509) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q16555
Expiration Date	12 months from date of receipt.