



Product Datasheet

MAP3K7 (Phospho-Thr187) Antibody [Out of stock] (orb191688)


Description	MAP3K7 (Phospho-Thr187) Antibody
Species/Host	Rabbit
Tested Applications	WB
Storage	Store at -20°C/1 year
Note	For research use only
Clonality	Polyclonal
Uniprot ID	O43318
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunohistochemical staining of human TB...



IHC-P staining of mouse testis tissue MA...