



Product Datasheet

Doublecortin antibody (orb19169)

Description	Goat polyclonal antibody to Doublecortin
Species/Host	Goat
Reactivity	Canine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	DCDC2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Full length recombinant human DCDC2 expressed in COS7 (Lane A), non-transfectant COS7 (lane B). Approx 52kDa band observed (calculated MW of 52.8kDa according to NP_057440.1). Data kindly provided by Wang and LoTurco, University of Connecticut, USA. Recommended concentration 0.1-0.3 µg/ml.
Protein Sequence	DGEKANKDAEQKED
Clonality	Polyclonal
MW	52.8
Entrez	51473
Dilution Range	ELISA: 1:4000, WB: 0.1-0.3ug/ml, IHC-P: 5-10µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of COS7 cell lysat...



Immunohistochemical staining of Human Ce...