



Product Datasheet

SNX9 antibody (orb19170)

Description	Goat polyclonal antibody to SNX9
Species/Host	Goat
Reactivity	Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Target	SNX9
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:1000.WB: Approx 60kDa band observed in human heart and human placenta lysates (calculated MW of 66.6kDa according to NP_057308.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	HKGAIEKVKESDK
Clonality	Polyclonal
MW	66.6
Entrez	51429
Dilution Range	ELISA: 1:4000, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human placenta ...



Immunofluorescence analysis of U2OS cell...