



Product Datasheet

CCKAR antibody (orb19173)

Description	Goat polyclonal antibody to CCKAR
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	CCKAR
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.IHC: In paraffin embedded Human Pancreas, shows cytoplasmic staining of some cells in the islets of Langerhans. Recommended concentration, 10ug/ml.WB: Preliminary experiments in Human Pancreas lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml).
Protein Sequence	FEASQKKSARKERKPST
Clonality	Polyclonal
MW	47.8
Entrez	886
Dilution Range	ELISA: 1:128000, IHC-P: 2-4µg/ml
Expiration Date	12 months from date of receipt.



orb19173 (10ug/ml)
staining of paraffin ...



Immunohistochemical
staining of Human
Pa...