

Product Datasheet

Prostaglandin E receptor 4 antibody (orb19187)

Description Goat polyclonal antibody to Prostaglandin E receptor 4

Species/Host Goat

Reactivity Canine, Human

Conjugation Unconjugated

Tested Applications ELISA, WB

Target EP4 prostaglandin receptor

Preservatives Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage Aliquot and store at -20°C. Minimize freezing and thawing.

Note For research use only

Application notes ELISA: Peptide ELISA: antibody detection limit dilution 1:16000. WB: Approx 55kDa band observed in Human Testes lysates (calculated MW of 53.1kDa according to NP_000949.1). Recommended concentration: 0.1-0.3 µg/ml.

Protein Sequence LEREVSKNPDLQA


Clonality Polyclonal

MW 53.1

Entrez [5734](#)

Dilution Range ELISA: 1:16000, WB: 0.1-0.3 µg/ml

Expiration Date 12 months from date of receipt.



orb19187
(0.1ug/ml)
staining of
Human Te...