



Product Datasheet

Histone H3K27me1 Antibody (orb191872)


Description	Histone H3K27me1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q16695
Expiration Date	12 months from date of receipt.



Dot blot analysis of all sorts of methyl...



Immunofluorescence analysis of 293T cell...



Western blot analysis of extracts of HeL...