

Product Datasheet

Histone H3K27me3 Antibody (orb191874)

Description	Histone H3K27me3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Form/Appearance	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1.0mg/ml
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Application notes	WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200
Clonality	Recombinant
MW	17KDa
Uniprot ID	Q16695
Expiration Date	12 months from date of receipt.