

Product Datasheet

CNBP antibody (orb19188)

Description	Goat polyclonal antibody to CNBP
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Target	ZNF9 / CNBP
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 19kDa band observed in Human Skeletal Muscle lysates (calculated MW of 19.5kDa according to NP_003409.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	GESGHLARECTIE
Clonality	Polyclonal
MW	19.5; 19.7; 19.6;19.0
Entrez	7555
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb19188
(2ug/ml)
staining of
Human
Musc...