

Product Datasheet

SORBS1 antibody (orb19190)

Description

Goat polyclonal antibody to SORBS1

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications
Target

SORBS1 / ponsin (internal)

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.IHC: In paraffin embedded Human Skeletal Muscle shows muscle fibres in longitudinal section. Recommended concentration, 2.5ug/ml.WB: Preliminary experiments in Human Adipose, Heart, Kidney, Liver and Muscle lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates

Protein Sequence

SPKSEDDDDSDLYSP

Clonality

Polyclonal

MW

76.6; 90.2; 142; 127; 101; 87.2; 83.0

Entrez


10580

Dilution Range

ELISA: 1:2000, WB: 1 µg/ml)., IHC-P: 5 µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of Human Sk...