


Product Datasheet

TIRAP antibody (orb19197)

Description	Goat polyclonal antibody to TIRAP
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TIRAP
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 24kDa band observed in lysates of Human myeloid erythroblast cell line K562 (calculated MW of 23.9kDa according to NP_001034750.1 and 24.8kDa according to NP_683708.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	QTLLKKPKRPNSE
Clonality	Polyclonal
MW	23.9; 24.8
Entrez	114609
Dilution Range	ELISA: 1:32000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb19197
(0.1µg/ml)
staining of
K562 cel...