



Product Datasheet

Interferon gamma antibody (orb19198)


Description	Goat polyclonal antibody to Interferon gamma
Species/Host	Goat
Reactivity	Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC
Target	AIF1 / IBA1 (isoform 3)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: Paraffin embedded Human Lung and Brain (Cortex). Recommended concentration: 5-6µg/ml.. Notes: Immunofluorescence: Strong expression of the protein seen in the cytoplasm and plasma membrane of U937, Caco-2 and U251 cells. Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of K562 cells
Protein Sequence	NKQFLDDPKYSSDED
Clonality	Polyclonal
MW	16.7
Entrez	199
Dilution Range	ELISA: 1:32000, WB: 0.5-1.5 µg/ml, FACS: 10ug/ml, IF/ICC: 10µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human frontal c...



Flow Cytometry analysis of K562 cells ti...



Immunofluorescence analysis of Caco-2 ce...