



Product Datasheet

CYP17A1 antibody (orb19207)

Description	Goat polyclonal antibody to CYP17A1
Species/Host	Goat
Reactivity	Canine, Feline, Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	CYP17A1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: In paraffin embedded Human Adenal Gland shows strong granular staining of the cytoplasm in medulla cells. Recommended concentration, 3-5µg/ml.WB: Approx 60kDa band observed in human ovary lysates (calculated MW of 57.4kDa according to NP_000093.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	LFKDGDKLEK
Clonality	Polyclonal
MW	57.4
Entrez	1586
Dilution Range	ELISA: 1:64000, WB: 0.1-0.3 µg/ml, IHC-P: 3-5 µg/ml
Expiration Date	12 months from date of receipt.



CYP17A1 antibody
(0.3 µg/ml) staining
o...



Immunohistochemical
analysis of formalin...