

Product Datasheet

VTI1B antibody (orb19209)

Description	Goat polyclonal antibody to VTI1B
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	VTI1B
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 26-28kDa band observed in Human Duodenum, Kidney and Testes lysates (calculated MW of 26.7kDa according to NP_006361.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	RDFDEKQQEANET
Clonality	Polyclonal
MW	26.7
Entrez	10490
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb19209
(1µg/ml)
staining of
Human
Duod...