

Product Datasheet

CYP3A4 antibody (orb19211)

Description	Goat polyclonal to CYP3A4
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC
Target	CYP3A4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: Paraffin embedded Human Liver. Recommended concentration: 5-7µg/ml.. Notes: Flow Cytometry: Flow cytometric analysis of HepG2 cells. Recommended concentration: 10ug/ml.
Protein Sequence	KESRLEDTQKHR
Clonality	Polyclonal
MW	57.3
Entrez	1576
Dilution Range	ELISA: 1:8000, IHC-P: 5-7 µg/ml
Expiration Date	12 months from date of receipt.