



Product Datasheet

GALR3 antibody (orb19212)

Description	Goat polyclonal antibody to GALR3
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Galanin Receptor 3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 37kDa band observed in Human Skin and Human Testes lysates (calculated MW of 39.6kDa according to NP_003605.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	QGPEPREGPVHGGE
Clonality	Polyclonal
MW	39.6
Entrez	8484
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml, IHC-P: 4-6µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Skin lysa...



Immunohistochemical staining of Human Te...