

Product Datasheet

WNV NS1(His Tag) (orb1921714)

Description	Flaviviruses belong to the family of flaviviridae, most of which are arboviruses mainly including Zi...
Tested Applications	CLIA, IA, LF
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Sample Types	Test sample type: Whole Blood. Test sample source: Human
Purity	> 90% pure (SDS-PAGE).
Source	HEK293
MW	39.61 kDa
Expiration Date	6 months from date of receipt.