

Product Datasheet

EBNA1 full length (orb1921728)

Description	Epstein-Barr virus (EBV) is a prevalent gamma herpesvirus, infecting upwards of 90% of the world pop...
Tested Applications	CLIA, IA, LF
Form/Appearance	Purified, Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Sample Types	Test sample type: Whole Blood. Test sample source: Human
Purity	> 90% pure (SDS-PAGE).
Source	HEK293
MW	56.43 kDa
Expiration Date	6 months from date of receipt.