

Product Datasheet

Mala P.f-mAb (Pairs) Antibody (orb1921740)

Description	Plasmodium is an obligate intracellular parasite belonging to the protozoan phylum apicomplexa. It s...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	LF
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Application notes	Recommended Applications: LFIA. Test Methods: Self-produced antigen test; Test 10 Positive Reference Materials; Test Sample Extraction Solution, Test 10 normal Nasal Swab Samples
Sample Types	Test sample type: Whole Blood. Test sample source: Human
Clonality	Monoclonal
Purity	> 95% pure (SDS-PAGE).
MW	1000 kDa
Expiration Date	12 months from date of receipt.