

Product Datasheet

S100 β Antigen (orb1921752)

Description	S100 β Antigen
Tested Applications	CLIA, LF
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Application notes	Recommended Applications: LFIA , CLIA. Test Methods: Use self-produced antigen test assay test 1ng/mL A05901; Use self-produced antigen test assay test 100 ng/ml A05901 (5 tests); Use our self-produced antigen test assay test 5 ng/mL A05901 (5 tests)
Sample Types	Test sample type: Serum / Plasma. Test sample source: Human
Purity	> 90% pure (SDS-PAGE).
Source	E. coli
Expiration Date	6 months from date of receipt.