

Product Datasheet

Soluble Growth Stimulation Expressed Gene2 (ST2) Antigen (orb1921755)

Description	Soluble Growth Stimulation Expressed Gene2 (ST2) Antigen
Tested Applications	CLIA, LF
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Application notes	Recommended Applications: LFIA , CLIA. Test Methods: Use self-produced antigen test assay test 1ng/mL A05801; Use self-produced antigen test assay test 100 ng/ml A05801 (5 tests); Use our self-produced antigen test assay test 5 ng/mL A05801 (5 tests)
Sample Types	Test sample type: Serum / Plasma. Test sample source: Human
Purity	> 90% pure (SDS-PAGE).
Source	HEK293
MW	36kDa
Expiration Date	6 months from date of receipt.