

Product Datasheet

CKMB-mAb (Pairs) Antibody (orb1921762)

Description	CKMB-mAb (Pairs) Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	CLIA, LF, Turbidimetry
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Application notes	Recommended Applications: LFIA, CLIA, TIA. Test Methods: Test Linear Range, Test Detection Limit
Sample Types	Test sample type: Serum / Plasma. Test sample source: Human
Clonality	Monoclonal
Purity	> 95% pure (SDS-PAGE).
MW	150 kDa
Expiration Date	12 months from date of receipt.