

## Product Datasheet

# Creatine Kinase Isoenzymes-MB (CKMB) Antigen (orb1921763)

|                            |                                                                        |
|----------------------------|------------------------------------------------------------------------|
| <b>Description</b>         | Creatine Kinase Isoenzymes-MB (CKMB) Antigen                           |
| <b>Tested Applications</b> | CLIA, LF, Turbidimetry                                                 |
| <b>Form/Appearance</b>     | Purified,Liquid                                                        |
| <b>Storage</b>             | For long term storage, the product should be stored at -20°C or lower. |
| <b>Note</b>                | For research use only                                                  |
| <b>Application notes</b>   | Recommended Applications: LFIA, CLIA, TIA.                             |
| <b>Sample Types</b>        | Test sample type: Serum / Plasma.<br>Test sample source: Human         |
| <b>Purity</b>              | > 90% pure (SDS-PAGE).                                                 |
| <b>Source</b>              | E. coli                                                                |
| <b>MW</b>                  | 85kDa                                                                  |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                         |