

Product Datasheet

Myo-mAb (Pairs) Antibody (orb1921765)

Description	Myoglobin (Mb) is a monomeric 17,800-Da protein, which incorporates a heme group in a hydrophobic cr...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	CLIA, LF
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Application notes	Recommended Applications: LFIA , CLIA. Test Methods: Test Linear Range, Test Detection Limit
Sample Types	Test sample type: Serum / Plasma. Test sample source: Human
Clonality	Monoclonal
Purity	> 95% pure (SDS-PAGE).
MW	150 kDa
Expiration Date	12 months from date of receipt.