

Product Datasheet

IL-6-mAb (Pairs) Antibody (orb1921772)

Description	Interleukin-6 (IL-6) is a 26 kDa secreted protein that consists of 184 amino acids with two N-glycos...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	CLIA, LF
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Application notes	Recommended Applications: LFIA , CLIA. Test Methods: Test Linear Range, Test Detection Limit
Sample Types	Test sample type: Serum / Plasma. Test sample source: Human
Clonality	Monoclonal
Purity	> 95% pure (SDS-PAGE).
MW	150 kDa
Expiration Date	12 months from date of receipt.