

Product Datasheet

IL6(His Tag) (orb1921774)

Description	Interleukin-6 (IL-6) is a 26 kDa secreted protein that consists of 184 amino acids with two N-glycos...
Tested Applications	CLIA, LF
Form/Appearance	Lyophilized from 0.22 µm filtered solution in PBS□5%mannital□0.01% Tween80 pH 7.4.
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Sample Types	Test sample type: Serum / Plasma. Test sample source: Human
Purity	≥ 95% as determined by reducing SDS-PAGE.
Source	CHO
MW	21.64 kDa
Expiration Date	6 months from date of receipt.