

Product Datasheet

SAA-mAb (Pairs) Antibody (orb1921775)

Description	Serum amyloid A (SAA) proteins comprise a family of molecules (SAA1-4), two members of which (SAA1 a...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	CLIA, LF, Turbidimetry
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Sample Types	Test sample type: Serum / Plasma. Test sample source: Human
Clonality	Monoclonal
Purity	> 95% pure (SDS-PAGE).
MW	150 kDa
Expiration Date	12 months from date of receipt.