

Product Datasheet

Transferrin-mAb (Pairs) Antibody (orb1921781)

Description	Transferrins (Tf) are glycoproteins found in vertebrates which bind to and consequently mediate the ...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	LF
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Sample Types	Test sample type: Stool. Test sample source: Human
Clonality	Monoclonal
Purity	> 95% pure (SDS-PAGE).
MW	150 kDa
Expiration Date	12 months from date of receipt.