

Product Datasheet

FOB-mAb (Pairs) Antibody (orb1921785)

Description	Fecal occult blood (FOB) test is clinically suitable for detecting bleeding in patients with gastroi...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	LF
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Application notes	Recommended Applications: LFIA. Test Methods: Self-produced antigen test; Test 10 Positive Reference Materials; Test Sample Extraction Solution, Test 10 normal Nasal Swab Samples; Test 500µg/mL Porcine Hemoglobin, Bovine Hemoglobin,Sheep Hemoglobin and Rabbit Hemoglobin
Sample Types	Test sample type: Stool. Test sample source: Human
Clonality	Monoclonal
Purity	> 95% pure (SDS-PAGE).
MW	150 kDa
Expiration Date	12 months from date of receipt.