

Product Datasheet

Helicobacter Pylori (H.pylori) Natural Antigen (orb1921792)

Description	Helicobacter pylori (H. pylori) is a gram-negative, microaerophilic, spiral-shaped and flagellated b...
Tested Applications	CLIA, LF
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Application notes	Recommended Applications: LFIA , CLIA. Test Methods: Use self-produced antigen test assay test 1ng/mL A99702; Use self-produced antigen test assay test 100 ng/ml A99702 (5 tests); Use our self-produced antigen test assay test 5 ng/mL A99702 (5 tests); Test Rotavirus,Adenovirus & Norovirus controls
Sample Types	Test sample type: Stool. Test sample source: Human
Purity	> 90% pure (SDS-PAGE).
Source	E. coli
Expiration Date	6 months from date of receipt.