

## Product Datasheet

### HPIV2-mAb (Pairs) Antibody (orb1921806)

|                            |                                                                                                                                                                                                                                                       |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | The human parainfluenza virus (HPIV) is an enveloped, negative-sense, single-stranded RNA virus that...                                                                                                                                               |
| <b>Species/Host</b>        | Mouse                                                                                                                                                                                                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                                          |
| <b>Tested Applications</b> | CLIA, LF                                                                                                                                                                                                                                              |
| <b>Form/Appearance</b>     | Purified,Liquid                                                                                                                                                                                                                                       |
| <b>Storage</b>             | For long term storage, the product should be stored at -20°C or lower.                                                                                                                                                                                |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                 |
| <b>Application notes</b>   | Recommended Applications: LFIA , CLIA.<br>Test Methods: Self-produced antigen test; Test 10 Positive Reference Materials; Test Sample Extraction Solution, Test 10 normal Nasal Swab Samples; Test Covid-19, Influenza A/B, Aendovirus & RSV controls |
| <b>Sample Types</b>        | Test sample type: Nasopharyngeal swab /Oropharyngeal swab / Nasal swab. Test sample source: Human                                                                                                                                                     |
| <b>Clonality</b>           | Monoclonal                                                                                                                                                                                                                                            |
| <b>Purity</b>              | > 95% pure (SDS-PAGE).                                                                                                                                                                                                                                |
| <b>MW</b>                  | 150 kDa                                                                                                                                                                                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                       |