

Product Datasheet

HPIV1-mAb (Pairs) Antibody (orb1921810)

Description	The human parainfluenza virus (HPIV) is an enveloped, negative-sense, single-stranded RNA virus that...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	CLIA, LF
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Sample Types	Test sample type: Nasopharyngeal swab /Oropharyngeal swab / Nasal swab. Test sample source: Human
Clonality	Monoclonal
Purity	> 95% pure (SDS-PAGE).
MW	150 kDa
Expiration Date	12 months from date of receipt.