

Product Datasheet

Flu B-mAb (Pairs) Antibody (orb1921815)

Description	Influenza B viruses are enveloped, negative-sense single-stranded viruses consisting of 8 genomic RN...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	CLIA, LF
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Application notes	Recommended Applications: LFIA , CLIA. Test Methods: Self-produced antigen test; Test 10 Positive Reference Materials; Test Sample Extraction Solution, Test 10 normal Nasal Swab Samples; Test Covid-19, Influenza A, Aendovirus & RSV controls
Sample Types	Test sample type: Nasopharyngeal swab /Oropharyngeal swab / Nasal swab. Test sample source: Human
Clonality	Monoclonal
Purity	> 95% pure(SDS-PAGE).
MW	1000 kDa
Expiration Date	12 months from date of receipt.