

Product Datasheet

Flu A-mAb (Pairs) Antibody (orb1921818)

Description	Influenza is an infectious respiratory disease; in humans, it is caused by influenza A (genus influe...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	CLIA, LF
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Application notes	Recommended Applications: LFIA , CLIA. Test Methods: Self-produced antigen test; Test 10 Positive Reference Materials; Test Sample Extraction Solution, Test 10 normal Nasal Swab Samples; Test Covid-19, Influenza B, Aendovirus & RSV controls
Sample Types	Test sample type: Nasopharyngeal swab /Oropharyngeal swab / Nasal swab. Test sample source: Human
Clonality	Monoclonal
Purity	> 95% pure(SDS-PAGE).
MW	150 kDa
Expiration Date	12 months from date of receipt.