

Product Datasheet

Recombinant Human FGF4 (orb1921836)

Description	Fibroblast Growth Factor 4 (FGF4) is one of a family of 18 secreted canonical FGF proteins that inte...
Endotoxins	≤10 EU/mg
Form/Appearance	Lyophilized Powder
Storage	Use a manual defrost freezer and avoid repeated freeze - thaw cycles. 12 months from date of receipt, -20 to -70°C as supplied. 1 month, 2 to 8°C under sterile conditions after reconstitution. 6 months, -20 to -87°C under sterile conditions after reconstitution.
Note	For research use only
Application notes	Application Notes: For research use only. Minimum Sample Size: 50µg/vial
Purity	≥95%
Source	CHO Cell
MW	19.2 kDa
Uniprot ID	P08620
Entrez	2249
Expiration Date	6 months from date of receipt.