

Product Datasheet

Recombinant Human IL-6 (orb1921848)

Description	Interleukin-6 (IL-6) is a 26 kDa secreted protein that consists of 184 amino acids with two N-glycos...
Endotoxins	≤10 EU/mg
Form/Appearance	Lyophilized Powder
Storage	Use a manual defrost freezer and avoid repeated freeze - thaw cycles. 12 months from date of receipt, -20 to -70°C as supplied. 1 month, 2 to 8°C under sterile conditions after reconstitution. 6 months, -20 to -74°C under sterile conditions after reconstitution.
Note	For research use only
Application notes	Application Notes: For research use only. Minimum Sample Size: 50µg/vial
Purity	≥95%
Source	CHO Cell
MW	20.64 kDa
Uniprot ID	P05231
Entrez	3569
Expiration Date	6 months from date of receipt.