

Product Datasheet

Transferrin - Antigen (orb1921853)

Description	Transferrins (Tf) are glycoproteins found in vertebrates which bind to and consequently mediate the ...
Endotoxins	≤0.5 EU/mg
Tested Applications	LF
Form/Appearance	Purified,Liquid
Storage	For long term storage, the product should be stored at -20°C or lower.
Note	For research use only
Application notes	Application Notes: For research use only. Minimum Sample Size: 100 mg/vial
Sample Types	Test sample type: Stool. Test sample source: Human
Purity	> 95% pure (SDS-PAGE).
Source	CHO
MW	75.2 kDa
Uniprot ID	P02787
Entrez	7018
Expiration Date	6 months from date of receipt.